

Patient Education Guide

Guidance and Support to
Help You Manage Your
Gastric-Jejunal Tube
(G-J Tube)



G-JET[®]

Low Profile Gastric-Jejunal
Enteral Feeding Tube



Table of Contents

G-JET[®]

Low Profile Gastric-Jejunal Enteral Feeding Tube

Applied Medical Technology, Inc. (AMT).....	3
Enteral Nutrition (EN).....	4
Introduction to Tube Feeding	6
How the G-JET Button is Different.....	8
G-JET Button Components.....	14
Attaching Feed & Extension Sets.....	16
Feeding Into the Jejunum.....	18
Decompressing the Stomach	20
Administering Medication	20
Dual Extension/Feed Set Hook-Up.....	21
Stoma Site Care & Concerns.....	22
Device Maintenance & Concerns.....	24
Feed/Extension Sets & Accessories.....	30
Available Sizes.....	32
AMT Clamp™ Device.....	33
Information Card.....	34



Applied Medical Technology, Inc. (AMT) is a global leader and manufacturer of enteral feeding devices and accessories, upholding the highest standard for performance and quality assurance. For over 25 years our core business has been enteral device development and manufacturing, engineering and designing our products in a state-of-the-art facility near Cleveland, Ohio. Our products include the MiniONE® line of low profile Balloon, Non-Balloon, and Capsule Non-Balloon buttons, the AMT G-JET® Low Profile Gastric-Jejunal Enteral Tube, and the AMT Bridle™ Nasal Tube Retaining System.


As a caregiver or patient, we believe you have a right to the most comfortable and reliable enteral tube. Referencing our Patient Education Guide will help provide you with the guidance and support to better understand enteral nutrition and our G-JET® device.

Enteral Nutrition

Proper nutrition is essential to maintaining our bodies' health, growth, and ability to heal. It is important to understand how the gastrointestinal (GI) or digestive system works in order to understand Enteral Nutrition (EN). The digestive system is made up of several hollow organs joined as a long tube or pathway through your body.

Digestion begins in the mouth, through chewing the food you intake into smaller pieces which are then swallowed and passed down the esophagus. The esophagus carries the food from your mouth to your stomach. Once in the stomach, the food gets broken down into its main nutrients. These nutrients then pass into the small intestine where they are absorbed into the blood stream. Once through the small intestine, the remaining material passes into the large intestine where excess liquid is absorbed leaving solid waste product behind to be passed out of the body.





Various medical conditions may make it difficult or impossible for a person to eat orally, denying the body of essential nutrients. In such cases, your physician will select the best type of alternative feeding for your situation.

A **gastrostomy tube**, also known as a g-tube, may be inserted through the abdominal wall to provide direct access to the stomach for feeding.

A **gastric-jejunal tube**, also known as a g-j tube, may be inserted through the abdominal wall into and through the stomach to provide direct access to the jejunum, which is the second part of the small intestine.

A g-tube or a g-j tube is a convenient and effective means for delivering nutritional formulas to the body. These nutritional formulas are either commercially available or homemade using a food processor. A physician will prescribe the proper feeding procedure, formula, and amount of water to most effectively feed you.

Introduction to Tube Feeding

A feeding tube, like the **AMT G-JET®** Button, is used for people who cannot take food through the mouth, or are unable to get the proper amount of nutrients to meet their current needs. The **AMT G-JET®** Button delivers a liquid food (called formula) directly into the small intestine by bypassing the mouth, the esophagus, and the stomach. Once the formula enters the body through the **AMT G-JET®** Button, digestion works the same as when food begins in the mouth. You get the same nutrients delivered to your body through a feeding tube as you would get by mouth.

Tube feeding with the **AMT G-JET®** Button provides you a complete and healthy diet so you can continue to grow and develop. It will take time to adjust to tube feeding at home, but with the help of your health care professionals, some practice, and this education guide, you will gain the confidence to take care of your **G-JET®** Button.

Additional Support:

In addition to your physician, there are many resources available to help you with Enteral Tube Feeding:

- **www.oley.org** - The Oley Foundation is for people who feed by tube or IV at home
- **www.parent-2-parent.com** - A unique support system for parents of children with developmental disabilities
- **www.feedingtubeawareness.com** - A group of parents of children with feeding tubes
- **www.foodfortubies.com** - An organization that helps people with feeding tubes use “real” food





Connect with AMT on Facebook®
<http://facebook.appliedmedical.net>




The AMT G-JET® Button

How it is different.

A Better Balloon:

Just like AMT's MiniONE® Balloon Buttons, the G-JET® Button's balloon is shaped like an apple to cover more surface area creating a better gastric seal to reduce leaks.



Device Meets biocompatibility standard ISO 10993 and is approved for permanent contact (greater than 30 days).



A Customizable Balloon:

Just like AMT's MiniONE® Balloon Buttons, the G-JET® Button's balloon fill volume can be easily adjusted for a custom-fit, reducing leaks and improving patient comfort.




People do not come in one size, so neither should their balloons! Your physician can recommend the best fill volume for you.

See page 24 for minimum, recommended, and maximum fill volumes.

A Durable Balloon:


- Balloon is made from the highest grade medical silicone
- More compliant and better able to withstand stomach contractions
- Device meets biocompatibility standard ISO 10993 (approved for permanent contact, greater than 30 days)



Reduces Risk of Misconnections:

Accidentally feeding the stomach instead of the jejunum can be dangerous! The G-JET® Button addresses this multiple ways.

A Clearly Labeled Device:



Jejunal length printed on "J" strap and laser engraved on the external bolster.

G and J Port Covers with safety plugs have distinctive shapes for tactile differentiation.




Glow Green Jejunal Access Feed Set




Mutually-Exclusive Ports:


- The Gastric Extension Set will only fit the gastric (white) port.
- The Jejunal Feed Set will only fit the jejunal (glow green) port.
- Device ports, extension and feed sets are color-coded.



Reduces Clogging:



The G-JET® Button transforms internally from a tri-lumen design in its proximal (gastric) segment to a single large lumen in its distal (jejunal) segment.




The G-JET's jejunal lumen provides a more than **60%** increase in cross sectional area to provide better flow.

With less clogging, patients may require fewer tube replacements.

Due to the larger J lumen, a guidewire is easy to remove. In addition, it is less likely to bend or kink the tube during guidewire removal.



Reduces Skin Irritation and Discomfort:



- 4 point stabilized external bolster creating more comfort and less risk of irritation
- No hard, bulky external bolster means more patient comfort




- Low profile and easier to hide under clothing
- Flexible bolster is easier to clean around and less likely to irritate the skin



Enhanced Features:

Anti-Kink Technology

Spring Reinforcement* maintains superior flow. Approved for MRI



Device is MRI Conditional



Easier Nighttime Connections

Jejunum port and jejunum connector Glow in the Dark

Improved Med Port

Higher plug retention reduces chances of “feeding the bed.”




*NOTE: Anti-kink spring only available in 16FR G-JET®

G-JET® Components

1 External Bolster

The external bolster maintains tube position. A properly fitted device should not fit tightly against the skin. Instead, an 1/8 inch (3mm) gap between the external base and the skin surface is recommended to provide room for air circulation and cleaning of the stoma site. The recommended amount of water to place into the balloon is printed on top of the fill valve.



2 Balloon Valve

The balloon is inflated and deflated by inserting a slip tip syringe into the balloon valve. **Do NOT attempt to feed through the balloon valve.** The recess in the valve can trap foreign material, such as lint or spilled formula. The balloon valve must be kept clean to function properly.

3 Anti-Reflux Valve / One Way Valve

A valve is located in the bottom of the feeding port (jejunal port) and the decompression/medication port (gastric port) to prevent stomach contents from leaking out of the tube. The valve is opened when the feeding set adapter is in place.

4 Internal Retention Balloon

A balloon is inflated inside the stomach to hold the tube in place. Once the tube is inserted, the balloon is filled with the prescribed amount of water. It is recommended to check the balloon volume as directed by your physician. See page 24 for complete instructions on how to check fill volume.

5 Feeding Port (Glow Green)

Used to administer formula into the jejunum. The feeding port should remain sealed by the port cover with safety plug when the device is not being used.




6 Decompression/Medication Port (White)

Used to decompress the stomach and administer medication as prescribed by your physician. The decompression/medication port should remain sealed by the port cover with safety plug when the device is not being used.

7 Individual Port Cover with Safety Plugs

Individual safety plugs for each port should be closed when not in use.



8 Lot Number & Device Size

The lot number and device size is etched (laser engraved) on the top of the device for complete traceability.

Attaching Feed & Extension Sets

Jejunal (Glow Green) Feed Set:

To Attach a Feed Set to the Jejunal Port of the G-JET® Button:

1. Remove Jejunal feeding port cover (see page 15, #7) from the top of the AMT G-JET® Button, revealing the glow green color interlock.
2. Holding the button, as shown in Figure A, line up the black mark on the Jejunal (glow green) Feeding Set Adapter with the black mark on the G-JET® device.
3. Make sure the feeding set adapter is pushed in completely and **gently** turn the adapter no more than 3/4 turn clockwise to lock the adapter in place. When fully locked, the adapter will stop rotating.




- **Do NOT attempt to rotate the adapter past 3/4 turn.**
If the adapter does not rotate freely when first inserted into the feeding port, make sure the adapter is **pushed in all the way** and that there is nothing blocking the jejunal port.
- Remove the feed set by rotating the connector COUNTER-CLOCKWISE until the black line of the feeding set aligns with the black line on the jejunal port. Remove the feeding set and cap the AMT Jejunal port with the attached feeding port cover.



Gastric (White) Extension Set:

To Attach an Extension Set to the Gastric Port of the G-JET® Button:

1. Remove the Gastric decompression/medication port cover (see page 15, #7) from the top of the AMT G-JET® Button, revealing the white color interlock.
2. Holding the button, as shown in Figure B, line up the black mark on the Gastric decompression/medication Extension Set with the black mark on the G-JET®.
3. Make sure the extension set adapter is pushed in completely and **gently** turn the adapter no more than 3/4 turn clockwise to lock the adapter in place. When fully locked, the adapter will stop rotating.



- **Do NOT attempt to rotate the adapter past 3/4 turn.**
If the adapter does not rotate freely when first inserted into the feeding port, make sure the adapter is **pushed in all the way** and that there is nothing blocking the gastric port.
- Remove the extension set by rotating the connector COUNTER-CLOCKWISE until the black line of the extension set aligns with the black line on the gastric port. Remove the set and cap the AMT Gastric port with the attached port cover.

Feeding into The Jejunum

For all feeding methods it is recommended that the patient be fed in an upright position, with their head up above the heart. The patient should be kept in an upright position for at least 30 minutes after feeding. Please consult your physician to create a feeding plan that meets your needs.

Jejunal Feeding:

- Gather the needed materials: Formula, Feed Set, Irrigating Syringe, Enteral Feeding Pump and Water for flushing the tube.
- **ALWAYS** wash your hands **BEFORE** preparing the formula and touching the G-JET[®] device. Wash with soap and warm water, rinse thoroughly and dry with a clean towel.
- Shake ready-to-use formula and wipe the container tops before opening with a clean cloth **OR** mix powdered or concentrated formula with water according to the nutritional label directions.
(ALWAYS follow the manufacturer's instructions. Throw away formula more than 24 hours old. **DO NOT** mix new formula with old formula.)
- **Prime the tubing:** Attach the feeding bag tubing to the right angle feed set and clamp one or both tubes. Pour the prescribed amount of formula into the feeding bag and allow the liquid to flow through to the end of the tubing where it is clamped off. This pushes the air out of the tubing and prevents it from being pumped into the stomach.
- Connect the Jejunal (Glow Green) Feed Set to the glow green port labeled "J" for jejunal. (For detailed direction see page 16)
- Connect feeding bag tubing to the pump. Set flow rate. Unclamp tubes and begin feeding. (Follow the manufacturer's directions to set the pump.)
- When feeding is complete, clamp both feed set tubing and feeding bag tubing. Disconnect feeding bag.
- Using Irrigating Syringe, flush feed set while it is connected to the G-JET[®] device with prescribed amount of water.
- Disconnect feed set and wash with warm soapy water.
- Never attempt to Bolus feed through the jejunal lumen.



Jejunal Feeding Continued:

- Flush the jejunal port EVERY SIX HOURS WITH AT LEAST 30 ml WATER. DO NOT USE FORCE. Note: If you are simultaneously decompressing the stomach while feeding and you see formula in the gastric drainage, stop the feeding and notify your specialist.
- **WARNING:** NEVER CONNECT THE JEJUNAL PORT TO SUCTION. DO NOT MEASURE RESIDUALS FROM THE JEJUNAL PORT.
- **CAUTION:** DO NOT ADMINISTER MEDICATIONS THROUGH THE JEJUNAL FEEDING PORT. THIS WILL CLOG THE TUBE. TO HELP AVOID CLOGGING, USE A FEEDING PUMP TO DELIVER FORMULA THROUGH THE JEJUNAL LUMEN. NEVER ALLOW FORMULA TO STAND IN THE TUBE.

Decompressing the Stomach

Your physician may recommend decompressing or “venting” the stomach (letting the air or liquid out) before or after feeding to help relieve pressure. ALWAYS FOLLOW YOUR SPECIALIST’S INSTRUCTIONS.

Stomach/Gastric Decompression:

- Connect the Gastric (White) Extension Set to the white port labeled “G” for gastric (For detailed directions see page 17). The end of the Extension Set may be left open for gravity drainage, or connected to low intermittent suction.
- Raise the Extension Set approximately 6" above the stomach so gas and stomach contents can pass freely up and down the tube. Do this as directed by your physician.
- DO NOT use continuous or high intermittent suction. High pressure could collapse the AMT G-JET® device or injure the stomach tissue and cause bleeding.
- Flush the gastric port EVERY SIX HOURS with at least 30 ml of water.



Administering Medication:

CAUTION: DO NOT administer medications through the jejunal feeding port.

- Use liquid medications when possible.
- Speak with your Pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, crush the solid or tablet medication into a fine powder using a pill crusher or the back of a spoon.
- Mix the powder with warm water until completely dissolved before administering through the Gastric (White) port of the AMT G-JET® device.
- DO NOT crush enteric-coated tablets, timed-release tablets or capsules.
- Give only ONE medication at a time and DO NOT mix medications.
- DO NOT add medications to food formulas or containers.
- Flush the gastric port EVERY SIX HOURS with at least 30 ml of water.

Dual Extension/Feed Set Hook-Up


Feed Set & Extension Set Hook-Up:

To Attach an Extension Set to the Gastric Port of the G-JET® Button:

- Follow the directions on page 17, steps 1 through 3.

To Attach a Feed Set to the Jejunal Port of the G-JET® Button:

- Follow the directions on page 16, steps 1 through 3. However, on step 3 only turn the feed set a 1/4 turn.



- **DO NOT attempt to rotate the adapter past 3/4 turn.** If the adapter does not rotate freely when first inserted into the gastric or jejunal port, make sure the adapter is **pushed in all the way** and that there is nothing blocking the gastric or jejunal port.
- Remove the extension set by rotating the connector COUNTER-CLOCKWISE until the black line of the extension or feed set aligns with the black line on the gastric or jejunal port. Remove the set and cap the AMT gastric or jejunal port with the attached port cover.

Stoma Site Care

- **The Stoma Site should be kept clean and dry at all times.** It is important to clean the stoma site daily with mild soap and warm water, or as directed by your physician.
- **DO NOT ROTATE THE TUBE.**
- Gather the needed materials to clean the stoma site: soap, clean cotton swabs or terry cloth, and warm water.
- **ALWAYS** wash your hands **BEFORE** touching the G-JET® device and the skin around the stoma site. Wash with soap and warm water, rinse thoroughly and dry with a clean towel.
- Clean the skin around the stoma site with clean cotton swabs or terry cloth using soap and warm water. Work from the device outwards in a circular motion.
- Clean the skin under the G-JET® device with a clean cotton swab or terry cloth using soap and warm water as directed by your physician.
- Rinse the whole site with warm water to remove all soap residue and allow the stoma site to completely air dry.
- **NOTE: BATH TIME.** Patients fitted with the G-JET® button are allowed to bathe and swim (make sure the safety plugs are in place). A good time for routine cleaning of the G-JET® device is during a bath. Make sure to allow the stoma site to completely air dry afterwards.





- **Always check the stoma site for redness, pain/soreness, swelling, or unusual drainage.** If any of these signs or symptoms are observed, contact a physician for advice.
- **Gauze or pads are *not* necessary.** If there is leakage, the G-J Tube may be too loose or too tight and should be remeasured. Contact your physician to have the stoma site remeasured.

Stoma Site Related Concerns:

Contact your Physician if any of the following is observed:

- If the skin around the stoma is red or raw
- If there is drainage around the site that is white, yellow, or green.
If this drainage smells bad
- If crusting is noted at the site
- If there is repetitive leakage of food or stomach contents (gauze or pads should not be necessary)
- If the patient experiences diarrhea or vomiting
- If the patient develops a fever
- Pain at the G-J stoma site
- Bleeding, puss or inflammation at the G-J stoma site
- If the tube falls out: ALWAYS have a spare G-TUBE, like AMT's MiniONE® Balloon Button, to place in the stoma tract that is similar in size to your G-JET® device. This will aid in keeping the stoma tract from closing. Once a replacement device is in place call your physician.

Device Maintenance

Balloon Maintenance:

Routine maintenance is recommended for the AMT G-JET® Button.

- Check the balloon fill volume as recommended by your physician.
- While holding the external bolster (see page 14, #1) in place, insert a syringe into the balloon fill valve (see page 14, #2) and withdraw the water and examine the amount.
- Refill the balloon with the originally prescribed amount of water or as prescribed by your physician.

Balloon Fill Volume Table:



French Size	Minimum	Recommended	Maximum
14FR	3 ml	4 ml	5 ml
16FR	3 ml	4 ml	5 ml



Feed Set & Extension Set Maintenance:

- Feed Set & Extension Sets and syringes should be cleansed with mild soap and warm water after each use.
- Be sure to rinse thoroughly to remove any formula or soap residual.
- For continuous feeding, the Feed Set should be cleaned at least once per day.
- To prevent clogging, always flush the Feed Set & Extension set after feeding and/or administering medications.



Contact your insurance provider regarding the number of Feed Set & Extension sets you are allowed per month.



Device Related Concerns:

Balloon Failure:

- Make sure the balloon is filled. (Resistance should be felt when gently pulling on the tube.) Make sure the balloon is filled to the prescribed volume, see the table on page 24.

NOTE: The recommended fill volume is printed on the balloon fill valve.

Balloon Will Not Deflate:

- Clean the balloon port with a clean cotton swab to make sure formula/medication or other contaminants are not blocking the balloon port (see page 14, #2). Insert a slip tip syringe, push and twist one-quarter turn. If the balloon still will not deflate, push the valve with a paper-clip to release the water. If the problem persists call your physician.

Constipation:

- May occur due to insufficient amounts of water in addition to formula, inactivity, change in formula or medication, or change in feeding routine. Call your doctor if this continues.

Diarrhea

- May occur if formula is spoiled or delivered too quickly. Mix new formula before each feeding. Deliver formula at a slower rate. Cleanliness is also very important. **All caregivers should wash their hands thoroughly prior to preparing the formula and handling the feeding sets.** Make sure the feeding sets are rinsed thoroughly to avoid soap in the formula. Call your doctor if diarrhea continues.



Feeding Tube Becomes Disconnected:

- Stop the feeding pump and estimate the amount of formula lost. Wipe tube connections thoroughly with soap and water. Clean the inside of the feed set and the device feeding port (see page 15, #5) with a clean cotton swab and soap and water. Dry connectors and reconnect tubes. Resume feeding, adding additional formula for the estimated loss.
- If disconnections persist, ask for an AMT Clamp™ device to eliminate disconnects. See page 31 for ordering details.

G-J Tube is Pulled Out of the Patient:

- Although the G-J Tube is designed to decrease the number of pull outs, a G-J Tube may become accidentally dislodged. ALWAYS have a spare G-TUBE, like AMT's MiniONE® Balloon Button, to place in the stoma tract that is similar in size to your G-JET® device. This will aid in keeping the stoma tract from closing. Once a replacement device is in place call your doctor.

G-J Tube is Too Tight Against the Patient's Skin:

- The G-J Tube should be able to move slightly from side to side easily without resistance from the patient's skin. If the G-J Tube does not move easily or redness/bleeding occurs at the stoma site or the area directly underneath the G-J Tube, call your doctor. The patient may need to be remeasured for a longer G-J Tube.

Leakage Around the Device:

- Make sure the balloon is filled. (Resistance should be felt when gently pulling on the tube.) Make sure the balloon is filled to the prescribed volume, see the table on page 24. ***NOTE: The recommended fill volume is printed on the top of the balloon fill valve.***
- Check that the G-J Tube corresponds to the prescribed French size and length. If the incorrect device is in place, call your doctor.
- The G-J Tube may be too tight or too loose. Call your doctor to have the stoma remeasured.
- For new placements, it may take time for the stoma tract to naturally heal, firm up around the tube, and conform to the balloon. If leaking persists contact your physician.

Causes of Leakage:

- Incorrect device size
- Weight change
- Amount of water in balloon may be too much or too little
- Tension at the stoma site (from extension sets)
- Stomach may need to be decompressed

Preventing Tube Blockage:

- Flush the G-J Tube with 5 ml (10-20 ml for adults) warm water before and after administering food or medication or every 6 hours for continuous feeding. ***DO NOT place foreign objects down the center of either port on the G-J Tube.*** This will damage the G-J Tube causing it to become nonfunctional.
- DO NOT mix medications with formula while feeding, always administer medication one at a time and thoroughly clean the extension set. Try to obtain your medications in liquid form when possible and if only available in tablet form ask your pharmacist if it is safe to crush to tablet into a fine powder and mix with warm water.



Upset Stomach:

- May occur if too much formula is administered or if the formula is delivered too quickly. Call you doctor if this continues.

Vomiting:

- Aspiration (inhalation of food or stomach contents into the lungs) may occur while vomiting, causing difficulty in breathing or other serious medical conditions. Stop feeding and decompress the stomach immediately if vomiting occurs.

Contact emergency care immediately if difficulty in breathing occurs or aspiration is suspected. DO NOT feed while the patient is laying flat.



Feed/Extension Set Accessories

G-JET® White Gastric Access Extension Set Options

Catalog No.	Length	Connector	Adapter
6-1221	12"	Right Angle	Bolus
6-2421	24"	Right Angle	Bolus
6-1211	12"	Straight	Bolus
6-2411	24"	Straight	Bolus
6-1222	12"	Right Angle	Y-Port (clear)
6-2422	24"	Right Angle	Y-Port (clear)
6-1212	12"	Straight	Y-Port (clear)
6-0211	2"	Straight	Bolus

G-JET® Glow Green Jejunal Access Feeding Set Options

Catalog No.	Length	Connector	Adapter
8-1257G	12"	Right Angle (Glow Green)	Y-Port (purple) Improved Med Port Plug
8-1255	12"	Right Angle	Y-Port (purple)
8-2455	24"	Right Angle	Y-Port (purple)
8-1221	12"	Right Angle	Bolus
8-2421	24"	Right Angle	Bolus
8-1222	12"	Right Angle	Y-Port (clear)
8-2422	24"	Right Angle	Y-Port (clear)
8-0211	2"	Straight	Bolus

For a complete list of products and accessories visit us on the web www.appliedmedical.net or call Customer Service at 800-869-7382



Connector and Adapter Options:



Gastric Right Angle Connector (white)



Jejunal Y-Port Adapter (purple)



Gastric Straight Connector (white)



Bolus Adapter



Jejunal Right Angle Connector (Glow Green)



Y-Port Adapter (clear)



Jejunal Straight Connector (beige)

Note: Feeding sets are supplied assembled. Choose the best combination of length, connector, and adapter type for your use.

Other Accessories

Catalog No.	Description
4-7000	Winged Luer Adapter (for direct G & J Port Access)
4-3000	AMT Clamp™ Device
1804	Guidewire Compatible Balloon Stoma Measuring Device



Winged Luer Adapter



Guidewire Compatible Balloon Stoma Measuring Device




AMT Clamp™ Device

Available G-JET Sizes

Cross Reference Guide:

AMT G-JET® Low Profile G-J Tube/
Mic-Key® Transgastric-Jejunal Feeding Tube



14 F		16 F	
AMT G-JET®	MIC-KEY®	AMT G-JET®	MIC-KEY®
GJ-1410-15	NA	GJ-1610-15	0270-16-1.0-15
GJ-1410-22	NA	GJ-1610-22	0270-16-1.0-22
GJ-1412-15	NA	GJ-1612-15	0270-16-1.2-15
GJ-1412-22	NA	GJ-1612-22	0270-16-1.2-22
GJ-1412-30	NA	GJ-1612-30	0270-16-1.2-30
GJ-1415-15	NA	GJ-1615-15	0270-16-1.5-15
GJ-1415-22	NA	GJ-1615-22	0270-16-1.5-22
GJ-1415-30	NA	GJ-1615-30	0270-16-1.5-30
GJ-1415-45	NA	GJ-1615-45	0270-16-1.5-45
GJ-1417-15	NA	GJ-1617-15	0270-16-1.7-15
GJ-1417-22	NA	GJ-1617-22	0270-16-1.7-22
GJ-1417-30	NA	GJ-1617-30	0270-16-1.7-30
GJ-1417-45	NA	GJ-1617-45	0270-16-1.7-45
GJ-1420-22	NA	GJ-1620-22	0270-16-2.0-22
GJ-1420-30	NA	GJ-1620-30	0270-16-2.0-30
GJ-1420-45	NA	GJ-1620-45	0270-16-2.0-45
GJ-1423-30	NA	GJ-1623-30	0270-16-2.3-30
GJ-1423-45	NA	GJ-1623-45	0270-16-2.3-45
GJ-1425-30	NA	GJ-1625-30	0270-16-2.5-30
GJ-1425-45	NA	GJ-1625-45	0270-16-2.5-45
GJ-1427-30	NA	GJ-1627-30	0270-16-2.7-30
GJ-1427-45	NA	GJ-1627-45	0270-16-2.7-45
GJ-1430-45	NA	GJ-1630-45	0270-16-3.0-45
GJ-1435-45	NA	GJ-1635-45	0270-16-3.5-45

How to Read Our Order Numbers

GJ - XXYY - ZZ

|
 French
 Size
 |
 Stoma
 Length
 |
 Jejunal
 Length


The AMT Clamp™ Device

Eliminate Pullouts with the AMT Clamp™ Device for Continuous Feeding

Our new Y-Port is made from a special elastomer that helps secure the stepped connector with a “sticky” friction fit.

The new purple color signifies “enteral feeding” to avoid misconnections, and the translucent material aids visualization so you can see the flow through the Y-Port walls.

If pullouts persist, ask for an AMT Clamp™ device to totally eliminate disconnects.



The AMT Clamp Device part number is 4-3000.

Information Card

Patient Information

Patient Name: _____

Physician Name: _____ Phone: _____

Placement Date: _____

Replacement Date: _____

G-JET® Device Information



Product No.: _____ Lot No.: _____

French Size: _____ Stoma Length (cm): _____

Jejunal Length: _____

Balloon Fill Volume (ml): _____

PLACE STICKER HERE

Fill in the above information or peel off die-cut sticker from tray.

Formula

Type of Formula: _____

Medications: _____

Homecare Company Name: _____

Phone: _____



Type of Feeding

Bolus: Y / N _____ Continuous: Y / N _____

Water to prime tube before feeding (ml): _____ Flow Rate (formula/hour): _____

Length of Feeding (minutes): _____ Additional Water (ml): _____

Additional Water (ml): _____

Special Instructions

Contact Information

Physician: _____ Phone: _____

HCP: _____ Phone: _____

**Save the MRI
Safety Card for
Your Records**



Patient Education Guide

G-JET® Gastric-Jejunal Feeding Tube

Additional Support:

In addition to your physician, there are many resources available to help you with Enteral Tube Feeding:

- **www.oley.org**
The Oley Foundation is for people who feed by tube or IV at home
- **www.parent-2-parent.com**
A unique support system for parents of children with developmental disabilities
- **www.feedingtubeawareness.com**
A group of parents of children with feeding tubes
- **www.foodfortubies.com**
An organization that helps people with feeding tubes use “real” food



A Better G-J Button:

- Less clogging
- Increased safety
- Improved comfort



The Leader in Enteral Device Innovation

An ISO 13485 Registered Company 

8006 Katherine Boulevard
Brecksville, OH 44141

Toll-Free: 800-869-7382

E-mail: cs@appliedmedical.net

www.appliedmedical.net

Patents Pending

Balloon Mini ONE® Button
US Patent numbers 6,019,746 and
D418,220 and additional patents pending

AMT Clamp™ Product
US Patent number 6,375,231

C-4418-Rev B

G-JET and Mini ONE are registered trademarks and AMT Clamp is a trademark of Applied Medical Technology, Inc. Mic-Key is a registered trademark of Kimberly-Clark Worldwide, Inc. © 2014 Applied Medical Technology, Inc.



Connect with AMT on Facebook®
<http://facebook.appliedmedical.net>



Made in USA